

## Capnography/Capnometry During Mechanical Ventilation: 2011

Brian K Walsh RRT-NPS FAARC, David N Crotwell RRT-NPS, and  
Ruben D Restrepo MD RRT FAARC

We searched the MEDLINE, CINAHL, and Cochrane Library databases for articles published between January 1990 and November 2010. The update of this clinical practice guideline is based on 234 clinical studies and systematic reviews, 19 review articles that investigated capnography/capnometry during mechanical ventilation, and the 2010 American Heart Association Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care. The following recommendations are made following the Grading of Recommendations Assessment, Development, and Evaluation (GRADE) scoring system: (1) Continuous-waveform capnography is recommended, in addition to clinical assessment to confirm and monitor correct placement of an endotracheal tube. (2) If waveform capnography is not available, a non-waveform exhaled CO<sub>2</sub> monitor, in addition to clinical assessment, is suggested as the initial method for confirming correct tube placement in a patient in cardiac arrest. (3) End-tidal CO<sub>2</sub> (P<sub>ETCO<sub>2</sub></sub>) is suggested to guide ventilator management. (4) Continuous capnometry during transport of the mechanically ventilated patients is suggested. (5) Capnography is suggested to identify abnormalities of exhaled air flow. (6) Volumetric capnography is suggested to assess CO<sub>2</sub> elimination and the ratio of dead-space volume to tidal volume (V<sub>D</sub>/V<sub>T</sub>) to optimize mechanical ventilation. (7) Quantitative waveform capnography is suggested in intubated patients to monitor cardiopulmonary quality, optimize chest compressions, and detect return of spontaneous circulation during chest compressions or when rhythm check reveals an organized rhythm. *Key words: capnography; capnometry; colorimetric CO<sub>2</sub>; end-tidal carbon dioxide; volumetric CO<sub>2</sub>.* [Respir Care 2011;56(4):503–509. © 2011 Daedalus Enterprises]

### CO<sub>2</sub> MV 1.0 DESCRIPTION

For the purposes of this clinical practice guideline, capnography refers to the evaluation of the CO<sub>2</sub> in the respiratory gases of mechanically ventilated patients. A capnographic device incorporates one of 2 types of sampling

---

Brian K Walsh RRT-NPS FAARC is affiliated with the Respiratory Care Department, Children's Medical Center, Dallas, Texas. David N Crotwell RRT-NPS is affiliated with the Department of Respiratory Care, Seattle Children's Hospital and Regional Medical Center, Seattle, Washington. Ruben D Restrepo MD RRT FAARC is affiliated with the Department of Respiratory Care, The University of Texas Health Sciences Center at San Antonio, San Antonio, Texas.

Correspondence: Ruben D Restrepo MD RRT FAARC, Department of Respiratory Care, The University of Texas Health Sciences Center at San Antonio, 7703 Floyd Curl Drive, MSC 6248, San Antonio TX 78229. E-mail: restrepo@uthscsa.edu.

DOI: 10.4187/respcare.01175

techniques: mainstream or sidestream.<sup>1</sup> Mainstream technique inserts a sampling window into the ventilator circuit for measurement of CO<sub>2</sub>, whereas a sidestream analyzer samples gas from the ventilator circuit, and the analysis occurs away from the ventilator circuit. Analyzers utilize infrared, mass or Raman spectra, or a photoacoustic spectra technology.<sup>1,2</sup> Flow measuring devices are utilized in volumetric capnographs. Colorimetric CO<sub>2</sub> detectors are a form of mainstream sampling, but are simplistic. The colorimetric CO<sub>2</sub> detector has a pH-sensitive chemical indicator that undergoes color change with each inspiration and expiration, thus reflecting the change in CO<sub>2</sub> concentration. These devices start at baseline color when minimal CO<sub>2</sub> is present and undergo gradual color change with increasing CO<sub>2</sub> concentration.<sup>3</sup>

### CO<sub>2</sub> MV 2.0 PROCEDURE

Capnography is the continuous analysis and recording of the CO<sub>2</sub> concentration in respiratory gas. Although the

terms capnography and capnometry are sometimes considered synonymous, capnometry means only the measurement of CO<sub>2</sub> in respiratory gas (ie, analysis alone), without a continuous written record or waveform. Capnographic waveforms may be time-based or volume-based.<sup>4</sup>

### CO<sub>2</sub> MV 3.0 SETTING

Capnography can be performed by trained healthcare personnel in any setting in which mechanically ventilated patients are found.

### CO<sub>2</sub> MV 4.0 INDICATIONS

There are 3 broad categories of indications for capnography/capnometry: verification of artificial airway placement; assessment of pulmonary circulation and respiratory status; and optimization of mechanical ventilation.

**4.1 Verification of Artificial Airway Placement.** Even when the endotracheal tube is seen to pass through the vocal cords and tube position is verified by chest expansion and auscultation during mechanical ventilation, providers should obtain additional confirmation of airway placement with waveform capnography or an exhaled CO<sub>2</sub> or esophageal detector device.<sup>5</sup>

**4.1.1 Exhaled CO<sub>2</sub> detectors,** including colorimetric and non-waveform, reliably detect intratracheal placement in patients whose cardiac output is not exceedingly low or who have not had prolonged circulatory failure. Their use in prolonged cardiac arrest merits further study.<sup>5,6</sup>

**4.1.1.1** When waveform capnography is not available, these methods can be used in addition to clinical assessment as the initial method for confirming correct tube placement in a patient in cardiac arrest.

**4.1.2** Capnography may be used as an adjunct to determine that tracheal rather than esophageal intubation has occurred.<sup>4,7,8</sup>

**4.1.3** All intubations must be confirmed by some form of P<sub>ETCO<sub>2</sub></sub> measurement.<sup>5,9</sup>

**4.1.4** Effective ventilation through a supraglottic airway device such as the laryngeal mask airway (LMA) should result in a capnograph waveform during cardiopulmonary resuscitation (CPR), and after return of spontaneous circulation.<sup>5</sup>

**4.1.5** When feasible, monitoring P<sub>ETCO<sub>2</sub></sub> during chest compressions is encouraged.<sup>5</sup>

**4.1.5.1** If the P<sub>ETCO<sub>2</sub></sub> is < 10 mm Hg during CPR, the clinician should attempt to improve the quality of compressions.

**4.1.5.2** An abrupt and sustained increase in P<sub>ETCO<sub>2</sub></sub> is a sensitive indicator of return of spontaneous circulation.

**4.1.6** P<sub>ETCO<sub>2</sub></sub> monitoring is one of the objective standards required for monitoring patients in transport, to ensure integrity of the airway.<sup>6,10,11</sup>

**4.1.6.1** Providers should observe a consistent capnographic waveform with ventilation to confirm and monitor endotracheal tube placement in the field, in the transport vehicle, on arrival at the hospital, and after any patient transfer, to reduce the risk of unrecognized tube misplacement or displacement.<sup>5,12</sup>

**4.1.7** Capnography can be used to detect inadvertent airway intubation during gastric tube insertion.<sup>13</sup>

**4.1.8** Life-threatening airway disasters and ventilator disconnection can be averted with continuous capnography.<sup>14-16</sup>

**4.2 Assessment of Pulmonary Circulation and Respiratory Status.** Capnography assists in:

**4.2.1** Determining changes in pulmonary circulation and respiratory status sooner than pulse oximetry. In patients without lung disease, substantial hypercarbia may present before pulse oximetry notifies the clinician of a change in ventilation.<sup>14,17-20</sup>

**4.2.2** Monitoring the adequacy of pulmonary, systemic, and coronary blood flow,<sup>20,21</sup> as well as estimation of the effective (non-shunted) pulmonary capillary blood flow by a partial rebreathing method.<sup>22-24</sup>

**4.2.3** Evaluating the partial pressure of exhaled CO<sub>2</sub>, especially P<sub>ETCO<sub>2</sub></sub>.

**4.2.4** Screening for pulmonary embolism.<sup>25-28</sup>

**4.3 Optimization of Mechanical Ventilation.** Capnography during mechanical ventilation allows:

**4.3.1** Continuous monitoring of the integrity of the ventilator circuit, including the artificial airway<sup>29</sup> or bag mask ventilation, in addition to potentially detecting mechanical ventilation malfunctions.<sup>30-32</sup>

**4.3.2** Decreasing the duration of ventilatory support.<sup>33</sup>

**4.3.3** Adjustment of the trigger sensitivity.<sup>34</sup>

**4.3.4** Evaluation of the efficiency of mechanical ventilation, by the difference between P<sub>aCO<sub>2</sub></sub> and the P<sub>ETCO<sub>2</sub></sub>.<sup>35</sup>

**4.3.5** Monitoring of the severity of pulmonary disease<sup>36,37</sup> and evaluating the response to therapy, especially therapies intended to improve the ratio of dead space to tidal volume (V<sub>D</sub>/V<sub>T</sub>) and ventilation-perfusion matching (V̇/Q̇).<sup>23,27,38-46</sup>

**4.3.6** Monitoring of V̇/Q̇ during independent lung ventilation.<sup>47,48</sup>

**4.3.7** Monitoring of inspired CO<sub>2</sub> when it is being therapeutically administered.<sup>49</sup>

**4.3.8** Graphic evaluation of the ventilator-patient interface. Evaluation of the capnogram may be useful in detecting rebreathing of CO<sub>2</sub>, obstructive pulmonary disease, the presence of inspiratory effort during neuromuscular blockade (curare cleft), cardiogenic oscillations, esophageal intubation, and cardiac arrest.<sup>50</sup>

**4.3.9** Measurement of the volume of CO<sub>2</sub> elimination to assess metabolic rate and/or alveolar ventilation.<sup>43,51-53</sup>

**4.3.10** Monitoring of V<sub>D</sub>/V<sub>T</sub> to determine eligibility for extubation in children.<sup>40,54</sup>

**4.3.11** There is a relationship between V<sub>D</sub>/V<sub>T</sub> and survival in patients with the acute respiratory distress syndrome.<sup>55-57</sup>

## CO<sub>2</sub> MV 5.0 CONTRAINDICATIONS

There are no absolute contraindications to capnography in mechanically ventilated patients, provided that the data obtained are evaluated with consideration given to the patient's clinical condition.

## O<sub>2</sub> MV 6.0 HAZARDS/COMPLICATIONS

Capnography with a clinically approved device is a safe, noninvasive test, associated with few hazards in most populations. Hazards/complications are different for the 2 types of capnographic device.

### 6.1 Mainstream

**6.1.1** Dead Space. Adapters inserted into the airway between the airway and the ventilator circuit should have a minimal amount of dead space. This effect is inversely proportional to the size of the patient being monitored.<sup>44,58</sup>

**6.1.2** The addition of the weight of a mainstream adapter can increase the risk of accidental extubation in neonates and small children.<sup>58</sup>

### 6.2 Sidestream

**6.2.1** The gas sampling rate from some sidestream analyzers may be high enough to cause auto-triggering when flow-triggering of mechanical breaths is used. This effect is also inversely proportional to the size of the patient.<sup>58</sup>

**6.2.2** The gas sampling rate can diminish delivered V<sub>T</sub> in neonates and small patients while using volume targeted or volume controlled ventilation modes.<sup>58</sup>

## CO<sub>2</sub> MV 7.0 LIMITATIONS OF PROCEDURE OR DEVICE

Capnography, when performed using a device calibrated and operated as recommended by the manufacturer, has

few limitations. It is important to note that although the capnograph provides valuable information about the efficiency of ventilation (as well as perfusion), it is not a replacement or substitute for assessing the P<sub>aCO<sub>2</sub></sub>.<sup>4,41,59-61</sup> The difference between P<sub>ETCO<sub>2</sub></sub> and P<sub>aCO<sub>2</sub></sub> increases as dead-space volume increases.<sup>62</sup> In fact, the difference between the P<sub>aCO<sub>2</sub></sub> and P<sub>ETCO<sub>2</sub></sub> varies in the same patient over time.<sup>43,63-65</sup> Alterations in breathing pattern and V<sub>T</sub> may introduce error into measurements designed to be made during stable, steady-state conditions.<sup>51,52,66</sup> Interpretation of results must take into account the stability of physiologic variables such as minute ventilation, V<sub>T</sub>, cardiac output,  $\dot{V}/\dot{Q}$ , and CO<sub>2</sub> body stores. Certain situations may affect the reliability of the capnogram. The extent to which the reliability is affected varies somewhat among types of devices. Limitations include:

**7.1** The composition of the respiratory gas may affect the capnogram (depending on the measurement technology incorporated).

**7.1.1** The infrared spectrum of CO<sub>2</sub> has some similarities to the spectra of both oxygen and nitrous oxide.<sup>50</sup> A high concentration of either oxygen or nitrous oxide, or both, may affect the capnogram, so a correction factor should be incorporated into the calibration of any capnograph used in such a setting.<sup>59</sup>

**7.1.2** The reporting algorithm of some devices (primarily mass spectrometers) assumes that the only gases present in the sample are those that the device is capable of measuring. When a gas that cannot be detected by the mass spectrometer (such as helium) is present, the reported CO<sub>2</sub> values are incorrectly elevated in proportion to the concentration of the gas present.<sup>4,67</sup>

**7.2** The breathing frequency may affect the capnograph. A high breathing frequency may exceed the capnograph's response capabilities. The presence of high airway resistance, respiratory rate, or inspiratory-to-expiratory ratio may decrease the accuracy of the measurement obtained from a sidestream capnograph, compared to a mainstream capnograph.<sup>68,69</sup> In addition, a breathing frequency > 10 breaths/min affects different capnographs differently.<sup>67</sup>

**7.3** Contamination of the monitor or sampling system by secretions or condensate, a sample tube of excessive length, too high a sampling rate, or obstruction of the sampling chamber can lead to unreliable results.

**7.4** Use of filters between the patient airway and the capnograph's sampling line may lead to artificially low P<sub>ETCO<sub>2</sub></sub> readings.<sup>31,70</sup>

**7.5** The sensitivity for confirmation of endotracheal intubation by color change could range from 67% to 72%.<sup>71</sup>

**7.6** Clinical conditions associated with false negative readings include:

**7.6.1** Low cardiac output may cause a false negative result when attempting to verify endotracheal tube position in the trachea.<sup>72</sup>

**7.6.2** During CPR a positive test confirms placement of the ETT within the airway, whereas a negative test indicates either esophageal intubation or airway intubation with poor or absent pulmonary blood flow and requires an alternate means of confirmation of tube position.<sup>73-75</sup>

**7.6.3** When the endotracheal tube is in the pharynx and when antacids and/or carbonated liquids are present in the stomach, a false negative reading may be present. However, the waveform does not continue during subsequent breaths.<sup>76</sup>

**7.6.4** Elimination and detection of CO<sub>2</sub> can be dramatically reduced in patients with severe airway obstruction and pulmonary edema.<sup>77</sup>

**7.7** Clinical conditions associated with false positive readings include:

**7.7.1** Colorimetric CO<sub>2</sub> detectors may give a false positive if contaminated with acidic or CO<sub>2</sub>-filled gastric content, intratracheal medications such as epinephrine, extreme humidity, or the presence of trichloroethylene or chloroform anesthetics. Most require at least 6 breaths before a decision can be made.<sup>6,78</sup>

**7.7.2** Detection of CO<sub>2</sub> in expired gas after esophageal intubation as a result of prior bystander mouth-to-mouth ventilation may result in a false positive reading.<sup>79</sup>

**7.7.3** A transient rise in P<sub>ETCO<sub>2</sub></sub> after sodium bicarbonate administration is expected, but should not be misinterpreted as an improvement in quality of CPR or a sign of return of spontaneous circulation.<sup>5</sup>

**7.8** Inaccurate measurement of expired CO<sub>2</sub> may be caused by leaks or other clinical circumstances preventing collection of expired gases,<sup>80</sup> including:

**7.8.1** Leaks in the ventilator circuit.

**7.8.2** Leaks around the tracheal tube cuff, an uncuffed tube, or the mask, including LMA.

**7.8.3** Bronchopleural fistula.

**7.8.4** Dialysis or extracorporeal life support.

## CO<sub>2</sub> MV 8.0 ASSESSMENT OF NEED

Capnography is considered a standard of care during general anesthesia. The American Society of Anesthesiologists has suggested that capnography be available for patients with acute ventilatory failure on mechanical ventilatory support.<sup>81</sup> The American College of Emergency Physicians recommends capnography as an adjunctive

method to ensure proper endotracheal tube position.<sup>75</sup> The 2010 American Heart Association Guidelines for Cardio-pulmonary Resuscitation and Emergency Cardiovascular Care recommend capnography to verify endotracheal tube placement in all age groups.<sup>6</sup> Assessment of the need to use capnography with a specific patient should be guided by the clinical situation. The patient's primary cause of respiratory failure and the severity of his or her condition should be considered.

## CO<sub>2</sub> MV 9.0 ASSESSMENT OF OUTCOME

Results should reflect the patient's condition and should validate the basis for ordering the monitoring. Documentation of results (along with all ventilatory and hemodynamic variables available), therapeutic interventions, and/or clinical decisions made based on the capnogram should be included in the patient's chart.

## CO<sub>2</sub> MV 10.0 RESOURCES

**10.1** Equipment: the capnograph and accessories (eg, airway adapter, sampling tube, depending on capnograph). The capnograph should be calibrated as recommended by the manufacturer.

**10.2** Personnel: licensed or credentialed respiratory therapists or individuals with similar credentials (eg, MD, RN) who have the necessary training and demonstrated skills to correctly calibrate and evaluate the capnograph, assess the patient and the patient-ventilator system, and the ability to exercise appropriate clinical judgment.

## CO<sub>2</sub> MV 11.0 MONITORING

**11.1** During capnography the following should be considered and monitored:

**11.1.1** Ventilatory variables: V<sub>T</sub>, respiratory rate, PEEP, ratio of inspiratory-to-expiratory time, peak airway pressure, and concentrations of respiratory gas mixture.<sup>3,38,44,72,82</sup>

**11.1.2** Hemodynamic variables: systemic and pulmonary blood pressure, cardiac output, shunt, and  $\dot{V}/\dot{Q}$  imbalances.<sup>23,41,66</sup>

## CO<sub>2</sub> MV 12.0 FREQUENCY

Capnography (or, at least, capnometry) should be available during endotracheal intubation.<sup>14,19,83</sup> Capnography is not indicated for every mechanically ventilated patient; however, when it is used, the measurement period should be long enough to allow determination of the P<sub>aCO<sub>2</sub></sub>-P<sub>ETCO<sub>2</sub></sub> difference, to note changes in the P<sub>aCO<sub>2</sub></sub>-P<sub>ETCO<sub>2</sub></sub> difference

as a result of therapy, and to allow interpretation of observed trends.

### CO<sub>2</sub> MV 13.0 INFECTION CONTROL

No specific precautions are necessary, although standard precautions (as described by the Centers for Disease Control and Prevention)<sup>27</sup> and precautions designed to limit the spread of tuberculosis<sup>80,84</sup> should always be implemented during patient care.

**13.1** Reusable mainstream sensors should be subjected to high-level disinfection between patients, according to the manufacturer's recommendations.

**13.2** The external surface of the monitor should be cleaned as needed, according to manufacturer's recommendations.

### CO<sub>2</sub> MV 14.0 RECOMMENDATIONS

The following recommendations are given based on the Grading of Recommendations Assessment, Development, and Evaluation (GRADE) scoring system<sup>85,86</sup>:

**14.1** Continuous waveform capnography is recommended in addition to clinical assessment as the most reliable method of confirming and monitoring correct placement of an endotracheal tube. (1A)

**14.2** If waveform capnography is not available, a non-waveform exhaled CO<sub>2</sub> monitor in addition to clinical assessment is suggested as the initial method for confirming correct tube placement in a patient in cardiac arrest. (2B)

**14.3** P<sub>ETCO<sub>2</sub></sub> is suggested as a method to guide ventilator management. (2B)

**14.4** Continuous capnometry during transport of a mechanically ventilated patient is suggested. (2B)

**14.5** Capnography is suggested to identify abnormalities of exhaled air flow. (2B)

**14.6** Volumetric capnography is suggested to assess CO<sub>2</sub> elimination and V<sub>D</sub>/V<sub>T</sub> to optimize mechanical ventilation. (2B)

**14.7** Quantitative waveform capnography is suggested in intubated patients to monitor CPR quality, optimize chest compressions, and detect return of spontaneous circulation during chest compressions or when rhythm check reveals an organized rhythm. (2C)

### CO<sub>2</sub> MV 15.0 CLINICAL PRACTICE GUIDELINE IDENTIFYING INFORMATION AND AVAILABILITY

#### 15.1 Adaptation

Original publication: *Respir Care* 1995;40(12):1321-1324.

#### 15.2 Guideline Developers

American Association for Respiratory Care Clinical Practice Guidelines Steering Committee.

Ruben D Restrepo MD RRT FAARC, Chair, Department of Respiratory Care, The University of Texas Health Sciences Center at San Antonio, San Antonio, Texas.

Brian K Walsh RRT-NPS FAARC, Children's Medical Center. Dallas, Texas.

David N Crotwell RRT-NPS, Seattle Children's Hospital. Seattle, Washington.

#### 15.3 Source(s) of Funding

None

#### 15.4 Financial Disclosures/Conflicts of Interest

Dr Restrepo is a consultant and researcher for Oridion, which manufactures capnographs. Mr Walsh and Mr Crotwell have disclosed no conflicts of interest.

### REFERENCES

- Block FE Jr, McDonald JS. Sidestream versus mainstream carbon dioxide analyzers. *J Clin Monit* 1992;8(2):139-141.
- O'Flaherty D. *Capnometry*. London: BMJ Publishing Group; 1994: 21-54.
- Garey DM, Ward R, Rich W, Heldt G, Leone T, Finer NN. Tidal volume threshold for colorimetric carbon dioxide detectors available for use in neonates. *Pediatrics* 2008;121(6):e1524-1527.
- Hess D, Branson RD. Noninvasive respiratory monitoring equipment. Philadelphia: Lippincott; 1994:184-216.
- Neumar RW, Otto CW, Link MS, Kronick SL, Shuster M, Callaway CW, et al. Part 8: adult advanced cardiovascular life support: 2010 American Heart Association guidelines for cardiopulmonary resuscitation and emergency cardiovascular care. *Circulation*. 2010;122(18 Suppl 3):S729-S767.
- Goldberg JS, Rawle PR, Zehnder JL, Sladen RN. Colorimetric end-tidal carbon dioxide monitoring for tracheal intubation. *Anesth Analg* 1990;70(2):191-194.
- Wenzel V, Voelckel WG, Krismer AC, Mayr VD, Strohmenger HU, Baubin MA, et al. [The new international guidelines for cardiopulmonary resuscitation: an analysis and comments on the most important changes]. *Anaesthesist* 2001;50(5):342-357. *Article in German*.
- Rudraraju P, Eisen LA. Confirmation of endotracheal tube position: a narrative review. *J Intensive Care Med* 2009;24(5):283-292.
- Field JM, Hazinski MF, Sayre MR, Chameides L, Schexnayder SM, Hemphill R, et al. Part 1: executive summary: 2010 American Heart Association guidelines for cardiopulmonary resuscitation and emergency cardiovascular care. *circulation*. 2010;122(18 Suppl 3):S640-S656.
- Braman SS, Dunn SM, Amico CA, Millman RP. Complications of intrahospital transport in critically ill patients. *Ann Intern Med* 1987; 107(4):469-473.
- Singh S, Allen WD Jr, Venkataraman ST, Bhende MS. Utility of a novel quantitative handheld microstreamcapnometer during transport of critically ill children. *Am J Emerg Med* 2006;24(3):302-307.
- Silvestri S, Ralls GA, Krauss B, Thundiyil J, Rothrock SG, Senn A, Carter E, Falk J. The effectiveness of out-of-hospital use of continuous end-tidal carbon dioxide monitoring on the rate of unrecognized misplaced intubation within a regional emergency medical services system. *Ann Emerg Med* 2005;45(5):497-503.

13. Howes DW, Shelley ES, Pickett W. Colorimetric carbon dioxide detector to determine accidental tracheal feeding tube placement. *Can J Anaesth* 2005;52(4):428-432.
14. Poirier MP, Gonzalez Del-Rey JA, McAneney CM, DiGiulio GA. Utility of monitoring capnography, pulse oximetry, and vital signs in the detection of airway mishaps: a hyperoxemic animal model. *Am J Emerg Med* 1998;16(4):350-352.
15. Ahrens T, Sona C. Capnography application in acute and critical care. *AACN Clin Issues* 2003;14(2):123-132.
16. Joint Commission. Preventing ventilator-related deaths and injuries. Sentinel Event Alert 2002, Issue 25, February 2002. [http://www.jointcommission.org/assets/1/18/sea\\_25.pdf](http://www.jointcommission.org/assets/1/18/sea_25.pdf). Accessed February 8, 2011.
17. Roberts WA, Maniscalco WM. A novel cause of error in capnographic confirmation of intubation in the neonatal intensive care unit. *Pediatrics* 1995;95(1):140-142.
18. Hall D, Goldstein A, Tynan E, Braunstein L. Profound hypercarbia late in the course of laparoscopic cholecystectomy: detection by continuous capnometry. *Anesthesiology* 1993;79(1):173-174.
19. Roberts WA, Maniscalco WM, Cohen AR, Litman RS, Chhibber A. The use of capnography for recognition of esophageal intubation in the neonatal intensive care unit. *Pediatr Pulmonol* 1995;19(5):262-268.
20. Shibutani K, Muraoka M, Shirasaki S, Kubal K, Sanchala VT, Gupte P. Do changes in end-tidal P<sub>CO2</sub> quantitatively reflect changes in cardiac output? *Anesth Analg* 1994;79(5):829-833.
21. Levine RL, Wayne MA, Miller CC. End-tidal carbon dioxide and outcome of out-of-hospital cardiac arrest. *N Engl J Med* 1997;337(5):301-306.
22. de Abreu MG, Geiger S, Winkler T, Ragaller M, Pfeiffer T, Leutheuser D, et al. Evaluation of a new device for noninvasive measurement of nonshunted pulmonary capillary blood flow in patients with acute lung injury. *Intensive Care Med* 2002;28(3):318-323.
23. de Abreu MG, Quintel M, Ragaller M, Albrecht DM. Partial carbon dioxide rebreathing: a reliable technique for noninvasive measurement of nonshunted pulmonary capillary blood flow. *Crit Care Med* 1997;25(4):675-683.
24. vanHeerden PV, Baker S, Lim SI, Weidman C, Bulsara M. Clinical evaluation of the non-invasive cardiac output (NICO) monitor in the intensive care unit. *Anaesth Intensive Care* 2000;28(4):427-430.
25. Rodger M, Wells PS. Diagnosis of pulmonary embolism. *Thromb Res* 2001;103(6):V225-238.
26. Rodger MA, Jones G, Rasuli P, Raymond F, Djunaedi H, Bredeson CN, et al. Steady-state end-tidal alveolar dead space fraction and D-dimer: bedside tests to exclude pulmonary embolism. *Chest* 2001;120(1):115-119.
27. Bolyard EA, Tablan OC, Williams WW, Pearson ML, Shapiro CN, Deitchmann SD. Guideline for infection control in health-care personnel, 1998. *Hospital Infection Control Practices Advisory Committee. Infect Control Hosp Epidemiol* 1998;19(6):407-463. Erratum in: *Infect Control Hosp Epidemiol* 1998;19(7):493.
28. Rumpf TH, Krizmaric M, Grmec S. Capnometry in suspected pulmonary embolism with positive D-dimer in the field. *Crit Care* 2009;13(6):R196.
29. Spahr-Schopfer IA, Bissonnette B, Hartley EJ. Capnometry and the paediatric laryngeal mask airway. *Can J Anaesth* 1993;40(11):1038-1043.
30. Muniz AE. False-negative capnographic reading caused by a malfunctioning bag-valve-mask device resulting in a pneumomediastinum. *Resuscitation* 2008;78(3):378-380.
31. Hardman JG, Mahajan RP, Curran J. The influence of breathing system filters on paediatric capnography. *Paediatr Anaesth* 1999;9(1):35-38.
32. Kumar AY, Bhavani-Shankar K, Moseley HS, Delph Y. Inspiratory valve malfunction in a circle system: pitfalls in capnography. *Can J Anaesth* 1992;39(9):997-999.
33. Cheifetz IM, Myers TR. Respiratory therapies in the critical care setting. Should every mechanically ventilated patient be monitored with capnography from intubation to extubation? *Respir Care* 2007;52(4):423-438; discussion 438-442.
34. Thompson JE, Jaffe MB. Capnographic waveforms in the mechanically ventilated patient. *Respir Care* 2005;50(1):100-108; discussion 108-109.
35. Kerr ME, Zempsky J, Sereika S, Orndoff P, Rudy EB. Relationship between arterial carbon dioxide and end-tidal carbon dioxide in mechanically ventilated adults with severe head trauma. *Crit Care Med* 1996;24(5):785-790.
36. Bedford NM, Hardman JG. Predicting patients' responses to changes in mechanical ventilation: a comparison between physicians and a physiological simulator. *Intensive Care Med* 1999;25(8):839-842.
37. Ghamra ZW, Arroliga AC. Volumetric capnography in acute respiratory distress syndrome: is the era of day-to-day monitoring finally here? *Respir Care* 2005;50(4):457-458.
38. Engoren M. Efficacy of capnometry in ventilatory management of cardiac patients. *J Cardiothorac Vasc Anesth* 1993;7(5):538-540.
39. Hardman JG, Aitkenhead AR. Estimation of alveolar deadspace fraction using arterial and end-tidal CO<sub>2</sub>: a factor analysis using a physiological simulation. *Anaesth Intensive Care* 1999;27(5):452-458.
40. Hubble CL, Gentile MA, Tripp DS, Craig DM, Meliones JN, Cheifetz IM. Dead-space to tidal volume ratio predicts successful extubation in infants and children. *Crit Care Med* 2000;28(6):2034-2040.
41. Jellinek H, Hiesmayr M, Simon P, Klepetko W, Haider W. Arterial to end-tidal CO<sub>2</sub> tension difference after bilateral lung transplantation. *Crit Care Med* 1993;21(7):1035-1040.
42. Kallet RH, Daniel BM, Garcia O, Matthay MA. Accuracy of physiologic dead-space measurements in patients with acute respiratory distress syndrome using volumetric capnography: comparison with the metabolic monitor method. *Respir Care* 2005;50(4):462-467.
43. Russell GB, Graybeal JM. Reliability of the arterial to end-tidal carbon dioxide gradient in mechanically ventilated patients with multisystem trauma. *J Trauma* 1994;36(3):317-322.
44. Szaflarski NL, Cohen NH. Use of capnography in critically ill adults. *Heart Lung* 1991;20(4):363-372.
45. Taskar V, John J, Larsson A, Wetterberg T, Jonson B. Dynamics of carbon dioxide elimination following ventilator resetting. *Chest* 1995;108(1):196-202.
46. McSwain SD, Hamel DS, Smith PB, Gentile MA, Srinivasan S, Meliones JN, et al. End-tidal and arterial carbon dioxide measurements correlate across all levels of physiologic dead space. *Respir Care* 2010;55(3):288-293.
47. Cinnella G, Dambrosio M, Brienza N, Giuliani R, Bruno F, Fiore T, et al. Independent lung ventilation in patients with unilateral pulmonary contusion. Monitoring with compliance and EtCO<sub>2</sub>. *Intensive Care Med* 2001;27(12):1860-1867.
48. Colman Y, Krauss B. Microstreamcapnography technology: a new approach to an old problem. *J Clin Monit Comput* 1999;15(6):403-409.
49. Fatigante L, Cartei F, Ducci F, Marini C, Prediletto R, Caciagli P, et al. Carbogen breathing in patients with glioblastomamultiforme submitted to radiotherapy. Assessment of gas exchange parameters. *Acta Oncol* 1994;33(7):807-811.
50. Bhavani-Shankar K, Moseley H, Kumar AY, Delph Y. Capnometry and anaesthesia. *Can J Anaesth* 1992;39(6):617-632.
51. Brandi LS, Bertolini R, Santini L, Cavani S. Effects of ventilator resetting on indirect calorimetry measurement in the critically ill surgical patient. *Crit Care Med* 1999;27(3):531-539.

52. Brandi LS, Santini L, Bertolini R, Malacarne P, Casagli S, Baraglia AM. Energy expenditure and severity of injury and illness indices in multiple trauma patients. *Crit Care Med* 1999;27(12):2684-2689.
53. Sullivan KJ, Kissoon N, Goodwin SR. End-tidal carbon dioxide monitoring in pediatric emergencies. *Pediatr Emerg Care* 2005;21(5):327-332; quiz 333-325.
54. Wratney AT, Cheifetz IM. Extubation criteria in infants and children. *Respir Care Clin N Am* 2006;12(3):469-481.
55. Nuckton TJ, Alonso JA, Kallet RH, Daniel BM, Pittet JF, Eisner MD, Matthay MA. Pulmonary dead-space fraction as a risk factor for death in the acute respiratory distress syndrome. *N Engl J Med* 2002;346(17):1281-1286.
56. Lucangelo U, Bernabè F, Vatua S, Degrassi G, Villagrà A, Fernandez R, et al. Prognostic value of different dead space indices in mechanically ventilated patients with acute lung injury and ARDS. *Chest* 2008;133(1):62-71.
57. Raurich JM, Vilar M, Colomar A, Ibáñez J, Ayestarán I, Pérez-Bárcena J, Llompert-Pou JA. Prognostic value of the pulmonary dead-space fraction during the early and intermediate phases of acute respiratory distress syndrome. *Respir Care* 2010;55(3):282-287.
58. Jacobus C. Noninvasive monitoring in neonatal and pediatric care. St. Louis: Elsevier; 2009:137-146.
59. Hess D. Capnometry. New York: McGraw-Hill; 1998:377-400.
60. Isert P. Control of carbon dioxide levels during neuroanaesthesia: current practice and an appraisal of our reliance upon capnography. *Anaesth Intensive Care* 1994;22(4):435-441.
61. Laffon M, Gouchet A, Sitbon P, Guicheteau V, Biyick E, Duchalais A, et al. Difference between arterial and end-tidal carbon dioxide pressures during laparoscopy in paediatric patients. *Can J Anaesth* 1998;45(6):561-563.
62. Russell GB, Graybeal JM. End-tidal carbon dioxide as an indicator of arterial carbon dioxide in neurointensivecare patients. *J NeurosurgAnesthesiol* 1992;4(4):245-249.
63. Russell GB, Graybeal JM. The arterial to end-tidal carbon dioxide difference in neurosurgical patients during craniotomy. *Anesth Analg* 1995;81(4):806-810.
64. Seguin P, Bleichner JP, Branger B, Guillou YM, Feuillu A, Malle-dant Y. [The measurement of end-tidal carbon dioxide (PETCO<sub>2</sub>) is not a significant parameter to monitor in patients with severe traumatic brain injury.] *Can J Anaesth* 2001;48(4):396-400. *Article in French.*
65. Grenier B, Verchere E, Mesli A, Dubreuil M, Siao D, Vandendriessche M, et al. Capnography monitoring during neurosurgery: reliability in relation to various intraoperative positions. *Anesth Analg* 1999;88(1):43-48.
66. Gamma de Abreu M, Melo MF, Giannella-Neto A. Pulmonary capillary blood flow by partial CO<sub>2</sub> rebreathing: importance of the regularity of the respiratory pattern. *Clin Physiol* 2000;20(5):388-398.
67. Graybeal J, M. Relative agreement between Raman and mass spectrometry for measuring end-tidal carbon dioxide. *Respir Care* 1994;39(3):190-194.
68. McEvedy BA, McLeod ME, Kirpalani H, Volgyesi GA, Lerman J. End-tidal carbon dioxide measurements in critically ill neonates: a comparison of side-stream and mainstream capnometers. *Can J Anaesth* 1990;37(3):322-326.
69. Tingay DG, Stewart MJ, Morley CJ. Monitoring of end tidal carbon dioxide and transcutaneous carbon dioxide during neonatal transport. *Arch Dis Child Fetal Neonatal Ed* 2005;90(6):F523-F526.
70. Hardman JG, Curran J, Mahajan RP. End-tidal carbon dioxide measurement and breathing system filters. *Anaesthesia* 1997;52(7):646-648.
71. Keller WR, Biehler J, Linares MY, Garcia-Pena BM. False-positive colorimetric capnometry after ingestion of carbonated beverages. *Pediatr Emerg Care* 2009;25(2):69-73.
72. Li J. Capnography alone is imperfect for endotracheal tube placement confirmation during emergency intubation. *J Emerg Med* 2001;20(3):223-229.
73. 2005 American Heart Association (AHA) guidelines for cardiopulmonary resuscitation (CPR) and emergency cardiovascular care (ECC) of pediatric and neonatal patients: pediatric advanced life support. *Pediatrics* 2006;117(5):e1005-e1028.
74. The International Liaison Committee on Resuscitation (ILCOR) consensus on science with treatment recommendations for pediatric and neonatal patients: pediatric basic and advanced life support. *Pediatrics* 2006;117(5):e955-e977.
75. American College of Emergency Physicians, Clinical Policies Committee. Verification of endotracheal tube placement. *Ann Emerg Med* 2002;40(5):551-552.
76. Sum Ping ST, Mehta MP, Symreng T. Accuracy of the FEF CO<sub>2</sub> detector in the assessment of endotracheal tube placement. *Anesth Analg* 1992;74:415-419.
77. Ward KR, Yealy DM. End-tidal carbon dioxide monitoring in emergency medicine. Part 2: clinical applications *Acad Emerg Med* 1998;5(6):637-646.
78. Kamlin CO, O'Donnell CP, Davis PG, Morley CJ. Colorimetric end-tidal carbon dioxide detectors in the delivery room: strengths and limitations. A case report. *J Pediatr* 2005;147(4):547-548.
79. Kramer-Johansen J, Dorph E, Steen PA. Detection of carbon dioxide in expired air after oesophageal intubation; the role of bystander mouth-to-mouth ventilation. *Acta Anaesthesiol Scand* 2008;52(1):155-157.
80. Kallet RH. Capnography and respiratory care in the 21st century. *Respir Care* 2008;53(7):860-861.
81. American Society of Anesthesiologists. Standards for basic anesthetic monitoring. <http://www.asahq.org/knowledge-base/ethics-and-medicolegal-issues/asa/~media/for%2520members/documents/standards%2520guidelines%2520stmts/basic%2520anesthetic%2520monitoring%25202011.ashx>. Accessed February 8, 2011.
82. Gentile MA, Cheifetz IM. Optimal positive end-expiratory pressure: the search for the Holy Grail continues. *Crit Care Med* 2004;32(12):2553-2554.
83. Sum Ping ST, Mehta MP, Symreng T. Reliability of capnography in identifying esophageal intubation with carbonated beverage or antacid in the stomach. *Anesth Analg* 1991;73(3):333-337.
84. Jensen PA, Lambert LA, Iademarco MF, Ridzon R. Guidelines for preventing the transmission of Mycobacterium tuberculosis in health-care settings, 2005. *MMWR Recomm Rep* 2005;54(RR-17):1-141.
85. Guyatt GH, Oxman AD, Vist GE, Kunz R, Falck-Ytter Y, Alonso-Coello P, et al. GRADE: an emerging consensus on rating quality of evidence and strength of recommendations. *BMJ* 2008;336(7650):924-926.
86. Jaeschke R, Guyatt GH, Dellinger P, Schunemann H, Levy MM, Kunz R, et al. Use of GRADE grid to reach decisions on clinical practice guidelines when consensus is elusive (abstract). *BMJ* 2008;337(Suppl):A744.